



Power Supplies

Recommendation of Power Supplies suitable for specific electrophoresis units

Unit	for 1-2 units	for 3-4 units
MINIeasy	Roth-MINI, EV243	Roth-MINI, EV243, EV231
PROfessional I		
MINI		
PROfessional II	Roth-STANDARD, EV243, EV231	Roth-STANDARD, EV231, EV265
MIDI 1	Roth-STANDARD, EV231, EV265	Roth-BLOT, EV265
PROfessional III		
MIDI PLUS		
MIDI 1 large	EV265, EV261	Roth-BLOT, EV265, EV261
PROfessional IV		
MIDI 2		Roth-BLOT, EV261
PROfessional V		
MAXI	Roth-STANDARD, EV231	Roth-BLOT, EV265, EV261
VARIA 1		
PROfessional Acetat		Roth-STANDARD, EV265
MULTI	EV215	EV233
MINI vertical	Roth-STANDARD, EV243, EV231	Roth-STANDARD, EV231, EV265
MAXI vertical	Roth-BLOT, EV243, EV261	Roth-BLOT, EV265, EV215
Sequencing	EV233	
IEF	EV215	
MINI Blotting	Roth-MINI, EV243, EV231	EV231, EV261
MAXI Blotting	Roth-BLOT, EV231, EV261	Roth-BLOT, EV261, EV202
Electroelution	Roth-MINI, Roth-STANDARD, EV243, EV261	Roth-STANDARD, EV261, EV265
<p><i>When connecting several units, a higher voltage must be applied in order to attain the same intensity of current and output and, subsequently, the same separating capacity.</i></p>		

Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5
 76185 Karlsruhe
 Postfach 100121
 76231 Karlsruhe
 Telefon: +49 (0) 721/ 5606-0
 Telefax: +49 (0) 721/ 5606-149
 E-Mail: info@carlroth.com
 Internet: www.carlroth.com

gk 08.2009